Changes in the operating manual

This supplementary sheet is a supplement to the operating manual for the BioSpectrometer basic and not a substitute. Therefore, also read the operating manual of the BioSpectrometer before using the BioSpectrometer for the first time. The current version of the operating manual is available on the internet at www.eppendorf.com

1.2.1 Hazard icons

Caution, general hazard point, observe the operating manual.

2.1 Main illustration

The location of the name plate deviates from the main illustration. It is located on the bottom of the device on the rear left.

2.3.5 Device self-check

The factory setting for the self-test interval has been adjusted to "weekly".

During the self-test, ensure that the following requirements are met:

- Observe the ambient conditions (chapter 11.2).
- Avoid temperature fluctuations (for example, caused by open windows).
- After the device has been transported from a cool to a warmer environment, condensation can form inside the device. Wait at least three hours before connecting the device to the power supply.
- Make sure that the cuvette shaft is clean.

4.4/4.5 Connecting a printer/PC

- ▶ Only connect devices that meet the safety requirements defined in IEC 60950-1.

5.3.2 Measurement procedure and 6.4.7 print & export

Method step print & export

In the **print & export** method step the [Sample] softkey is active in the display.

- If you want to select sample results for export or printout, press the [Sample] softkey.

6.2.2 Routine method group

Proteins direct UV: 2 additional methods

<i>Protein A 280</i>	Absorbance values are displayed directly as a result.
Evaluation using factor = 1	

<i>Protein A 280 1 mm</i>	For microliter cuvettes with a light path of 1 mm, e.g., Hellma® TrayCell or microliter cuvettes of a similar design.
Parameter Cuvette = 1 mm	

Nucleic acids: 1 additional method

<i>dsDNA 1 mm</i>	For microliter cuvettes with a light path of 1 mm, e.g., Hellma® TrayCell or microliter cuvettes of a similar design.
Parameter Cuvette = 1 mm	

11.1 Technical data: power supply

Power supply	100 V to 240 V \pm 10 %, 50 Hz to 60 Hz
Overvoltage category	II
Degree of contamination	2
Power consumption	Maximum power consumption according to name plate: 25 W Approx. 15 W during operation Approx. 5 W with dimmed display
Protection class	I
Fuses	T 2.5 A/250 V, 5 mm \times 20 mm (2 pieces)

11.2 Technical data: ambient conditions

Atmospheric pressure	Use up to an altitude of 2000 m above MSL
----------------------	---

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Eppendorf AG · 22331 Hamburg · Germany · Tel: +49 40 53801-0 · Fax: +49 40 538 01-556 · E-mail: eppendorf@eppendorf.com
Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA