## Changes in the operating manual

This supplementary sheet is a supplement to the operating manual for the BioSpectrometer kinetic and not a substitute. Therefore, also read the operating manual of the BioSpectrometer before using the BioSpectrometer for the first time. The current version of the operating manual is available on the internet at [www.eppendorf.com](http://www.eppendorf.com)

### 1.2.1 Hazard icons

Caution, general hazard point, observe the operating manual.

### 2.1 Main illustration

The location of the name plate deviates from the main illustration. It is located on the bottom of the device on the rear left.

### 2.3.5 Device self-check

The factory setting for the self-test interval has been adjusted to "weekly".

During the self-test, ensure that the following requirements are met:

- Observe the ambient conditions (chapter 11.2).
- Avoid temperature fluctuations (for example, caused by open windows).
- After the device has been transported from a cool to a warmer environment, condensation can form inside the device. Wait at least three hours before connecting the device to the power supply.
- Make sure that the cuvette shaft is clean.

### 4.4/4.5 Connecting a printer/PC

- ▶ Only connect devices that meet the safety requirements defined in IEC 60950-1.

### 5.3.2 Measurement procedure and 6.4.7 print & export

#### Method step print & export

In the **print & export** method step the [Sample] softkey is active in the display.

- If you want to select sample results for export or printout, press the [Sample] softkey.

#### 6.2.1 Absorbance method group

**Single  $\lambda$  – continuous:** The measured values are displayed in an absorbance-time graph.

- To subsequently evaluate measured values as linear regression kinetics, press the [Next >] softkey and go to the **process results** method step.

#### 6.2.2 Routine method group

**Proteins direct UV:** 2 additional methods

<i>Protein A 280</i> Evaluation using factor = 1	Absorbance values are displayed directly as a result.
<i>Protein A 280 1 mm</i> Parameter <b>Cuvette</b> = 1 mm	For microliter cuvettes with a light path of 1 mm, e.g., Hellma® TrayCell or microliter cuvettes of a similar design.

**Nucleic acids:** 1 additional method

<i>dsDNA 1 mm</i> Parameter <b>Cuvette</b> = 1 mm	For microliter cuvettes with a light path of 1 mm, e.g., Hellma® TrayCell or microliter cuvettes of a similar design.
--	---

### 11.1 Technical data: power supply

Power supply	100 V to 240 V $\pm$ 10 %, 50 Hz to 60 Hz
Overvoltage category	II
Degree of contamination	2
Power consumption	Maximum power consumption according to name plate: 50 W Approx. 30 W during operation Approx. 5 W with dimmed display and temperature control switched off
Protection class	I
Fuses	T 2.5 A/250 V, 5 mm $\times$ 20 mm (2 pieces)

### 11.2 Technical data: ambient conditions

Atmospheric pressure	Use up to an altitude of 2000 m above MSL
----------------------	---

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Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA