

115161 Histosec® pastilles (without DMSO)

solidification point 56-58°C embedding agent for histology

For general questions please contact our
Customer Service:

Merck KGaA
Frankfurter Str. 250
64293 Darmstadt
Germany
Phone: +49 6151 72-0
Fax: +49 6151 72 2000

04 March 2014



Histosec® is an [embedding media](#) for histology. It is a [paraffin](#) enriched with polymers, with or without the addition of [DMSO](#) (dimethyl sulfoxide). DMSO increases the rate of penetration of paraffin and provides additional preservation. Its melting range is between 56 and 58 °C and, within this range, any specimen can be processed without fear of damage irrespective of the method used.

Product number	Packaging	Qty/Pk
1151612504	Plastic bag	4 x 2.5 kg
1151619025	Fibre carton	25 kg

Alternatives

[111609](#) [Histosec® pastilles solidification point 56-58°C embedding agent for histology](#)

Product information

HS Code | 2712 20 90

Chemical and physical data

Solubility	(20 °C) practically insoluble
Melting point	60 °C
Bulk density	400 kg/m ³
Boiling point	>350 °C
Flash point	>100 °C

Safety information

Storage class	10 - 13 Other liquids and solids
WGK	NWG not water endangering
Disposal	3 Relatively unreactive organic reagents should be collected in container A. If halogenated, they should be collected in container B. For solid residues use container C.